

How to Optimize Extract Usage with



For optimized usage, follow the steps below*:

- Fill an unused GREER® Pick® Well™ with 0.125 mL of extract
- Handle the GREER® Pick® by the larger end and immediately insert it into the well
- Remove the GREER Pick from the well and briefly inspect to ensure that a small droplet of extract is present at the tip of the device
- Hold the GREER Pick at an approximate 45-degree angle to the skin at the test site, bringing the tip in contact with the skin in the direction the GREER Pick is pointing
- Press the tines, making contact with the superficial skin
- Withdraw the GREER Pick with a slight lifting of the skin
- Discard the used device into a suitable biohazard container
- Repeat steps until testing is complete



*GREER® Pick® Directions for Use.