

From your insights came this innovation.

The Greer Versa Vial Rack™

“I’ve got a dozen different racks, with different vial sizes, all cluttering up my work area.”

Three Rack Styles. Multiple Vial Sizes Per Rack.

“Have you ever tried lifting a full rack off a counter? Without spilling?”

Four Handles. One On Each Side.

“We don’t have a lot of room in the refrigerator to waste.”


Small Footprint With Vertical Sides.
Stackable, Too.



Devoted to Allergy Immunotherapy

The Greer Versa Vial Rack™

What Customers Want	What The Versa Vial Rack Delivers
Specific types for specific uses	3 styles to choose from: <ul style="list-style-type: none"> • Prescription • Testing • Bulk
One rack that can hold many vial sizes	Tiered wells hold multiple vial sizes
Easy to move and hold even when full	<ul style="list-style-type: none"> • 4 easy-grip handles • Vials fit snugly in tiered wells
Easy to clean	<ul style="list-style-type: none"> • Drains in every well • Matte finish cleans with soapy water or 70% isopropyl alcohol
Durable	High-impact polystyrene
Space saving	2 different base sizes designed to fit most refrigerators
Resilient	Matte-finish plastic designed to resist scuffs and hide scratches
Neat	<ul style="list-style-type: none"> • Stackable • Consistent color system



Cross-section showing multisized tiers

Prescription Versa Vial Rack



	Item # VRP003
Size	8" x 10"
Vial Capacity	48
Vial Sizes	5-20 vials 10-13 vials 10-20 vials

VRP003

Testing Versa Vial Rack

	Item # VRT002
Size	8" x 10"
Vial Capacity	48
Vial Sizes	5 mL scratch vials 5-13 vials Labeled negative and positive controls in 2 sizes



VRT002

Bulk Versa Vial Rack



	Item # VRB001
Size	9 1/8" x 13 1/8"
Vial Capacity	39
Vial Sizes	50-20 vials 30-20 vials 10-13 vials

VRB001



Devoted to Allergy Immunotherapy

The Greer Versa Vial Rack™

Versatile

- 3 styles
 - Bulk
 - Testing
 - Prescription
- Tiered design to hold multiple vial sizes
- Flat, matte-finish sides for easy labeling and less smearing

Portable

- Easy-grip handles on all 4 sides
- Snug fit for vials

Space saving

- 2 base sizes to fit in most refrigerators
- Stackable when empty or filled with same-size, sealed vials*

Easy cleaning

- Drain holes in every well
- Matte finish

Durable

- High-impact, chemical-resistant polystyrene
- Scratch- and smudge-hiding matte finish

*Excluding testing vial rack



For more information,
please contact your
Greer Customer Care
Specialist at 800.378.3906
or visit www.greerlabs.com.